

ZUMMO 

# JUICE EXTRACTION KITS

---

www.zummocorp.com



## New colour kits

**For all Nature models.**

Zummo innovates in its juice extraction kits, providing them with different colours depending on their size. The objective is for them to be easier to identify, thus improving usability for the machine operators, and obtaining maximum yield from each fruit.

# 1.

## How do we identify the kits?

Each size is represented by a colour and is recommended for different fruit, depending on the size:

- **Kit S (Grey)**: for limes, tangerines and small oranges.
- **Kit M (Orange)**: for standard oranges and lemons.
- **Kit L (Maroon)**: for large oranges, grapefruit, and pomegranates.
- **Kit XL (Pink)**: for extra-large oranges, grapefruit, and pomegranates.



\*The kits in the picture correspond to the Z40 Nature model.

# 2.

## How to choose the right juice extraction kit to maximise yield?

You must always select the right kit for the fruit calibre. The fruit diameter may be compatible with two juice extraction kits. In that case, it is advisable to choose the smallest kit of the two.

For example, if we are going to squeeze 74 mm Ø oranges in Z14, we could either use kit L or kit M. In this case, we would recommend using kit M.

# 3.

## Comparison table

		S		M		L		XL	
Ø mm		40-60	45-60	53-76	55-75	64-88	65-90	70-90	77-100
Model	Z1	—	—	—	Standard	—	Optional	—	—
	Z06	—	Optional	—	Standard	—	—	Standard	—
	Z14	—	—	—	Standard	—	—	Standard	—
	Z40	Optional	—	Standard	—	Optional	—	—	Optional

## \*Notes

- Kits of the same colour share the size name (S, M, L or XL), but not the diameter or the juice extraction range, which may vary depending on the juicer model.
- The colour kits are not exchangeable between the different models. But they are exchangeable between the different variants of each juicer model.
- The colours of the kits are independent of the colours of the machine containers.